

SET-UP

- Yes
- No
- Regular
- Immediate
- Partial
- Articulated
 - Denar
 - Hanau
- Sam
- Whip-Mix
- Panadent
- Wafer Bite™

Sculpturing Instructions

- Remove all attachments
- Retain first molar bands
- Retain lower 3-3 retainer
- Remove lower 3-3 retainer



Reset

- All teeth
 - Teeth circled
- R

8	7	6	5	4	3	2	1
1	2	3	4	5	6	7	8

 L

Overbite

- Ideal (1-2 mm)
- Other

Overjet

- Ideal (0 mm)
- Other

Anterior Root Torque

- Upper**
 - Same
 - 2° Palatally
 - 2° Labially
- Lower**
 - Same
 - 2° Lingually
 - 2° Labially

Spaces

- Close all
- Compromise
- Leave space between

In case of discrepancy between upper and lower arches, I prefer

- Good Cl. I molar relation
- Good Cl. I cuspid relation
- Space between cuspid & bicuspid
- Space between lateral & cuspid

Occlusal Plane

- Flat
- Curve of spee

Arch Form

- Ideal (standard)
- Approx. same
- Straight arch

Arch Width

- Upper**
 - Same
 - Expand ___mm
 - Constrict ___mm
- Lower**
 - Same
 - Expand ___mm
 - Constrict ___mm

Models Enclosed

- Recent, appliances off
- Original
- Appliances on
- Completed set-up

Impressions Enclosed

- Appliances off
- Appliances on

CUSTOM PERFECTOR®

U.S. PAT. NO. 5,415,542 FOREIGN PATS. PENDING



Made from translucent silicone.


Labial Wire with Acrylic Shield for Maximum Control and Retention

- ### Labial Wire
- Cuspid to Cuspid
 - Bicuspid to Bicuspid
 - Other—indicate on drawing



Seating Springs

- Yes
 - No
- Use arrows to indicate locations



End Appliance Distal To

- Appliance should cover all teeth to prevent super-eruption.
- First molars
 - Second molars
 - Other _____

Hinge Axis

- Use average
- Headplate tracing enclosed
- Headplate enclosed
- Articulated Set-up — see left column

Note: Hinge-axis opening during construction should equal or surpass that of patient's freeway space. This maximizes patient comfort, helps insure continuous seating of lower arch and encourages exercise wearing.

Airways

- No
- Yes (Freeway space must be 2-3 mm)

Lower Anterior Relief

- Trim to allow for future use of lower Spring Aligner while wearing Perfector.
- Wax lingual of lower anteriors to allow for bonded 3-3 retainer.

Special Instructions _____

Centric Occlusion By

- Backs parallel
- Wax bite
- Lines on buccal surfaces of molars

Return Case By

- 1st class
- UPS (where available)
- Overnight

--	--	--	--	--	--	--	--



TP Orthodontics, Inc.

www.tportho.com

HEADQUARTERS

Mailing & Shipping
 100 Center Plaza
 La Porte, Indiana 46350-9672
 Phone: 800-348-8856
 219-785-2591

WESTERN UNITED STATES

Mailing
 P.O. Box 742
 Lodi, California 95241-0742

Shipping

130 N. Houston Lane
 Lodi, California 95240-2405
 Phone: 800-343-5997
 209-368-7545

ACCOUNT NO. _____

Dr. _____

Address _____

City _____ State _____

Zip _____ Phone _____

Patient's Name _____

Date Shipped To TP	Date Required
<input type="checkbox"/> Please note my address change	<input type="checkbox"/> This is my first case with TP.

PLEASE DO NOT WRITE IN THIS SPACE			
500-100 <input type="checkbox"/>	500-106 <input type="checkbox"/>	400-600 <input type="checkbox"/>	400-800 <input type="checkbox"/>
408-610 <input type="checkbox"/>	508-700 <input type="checkbox"/>	400-700 <input type="checkbox"/>	400-850 <input type="checkbox"/>
500-104 <input type="checkbox"/>	508-800 <input type="checkbox"/>	400-900 <input type="checkbox"/>	
500-102 <input type="checkbox"/>	509-000 <input type="checkbox"/>	450-035 <input type="checkbox"/>	

PLEASE SEND ADDITIONAL SUPPLIES

(Fill in address label only if additional material requested)

Dr. _____

Address _____

City _____

State _____ Zip _____

- Spring Aligner/Retainer Rx
- Perfector®/Set-up Rx
- Positioner/Set-up Rx
- Bondable Retainer & Indirect Bonding Rx
- Model Sculpture Rx
- Shipping Boxes
- Shipping Labels
- Shipping Bags
- Hinge-Axis Analysis Sheets
- Other

Return white and yellow copy. Retain pink copy.